

Robohand Information and Care

Please remove Robohand and inspect skin every hour for the first week of use. Inspect skin for chaffing, redness or any complaints of pain or numbness. If any of these conditions develop, please stop use and contact Robohand USA at (888)445-2178.

Please do not allow Custom Orthoplastic to get hot, it will lose its shape. Do not leave in hot car, place near stove, heater, windows or clothes dryers.

Robohand can be cleaned with cool water and soap, wiped with baby wipes or alcohol.

Ink can be removed with chlorine bleach on a cotton swab.

Velcro straps can be washed by hand or in washing machine (in sweater bag). Velcro straps and tabs will need replaced with time.

Please be sure both Robohand and wearers arm are very dry before using Robohand.

The Robohand can be used for part of the day initially and increased as tolerated to avoid overuse. Use can be increased gradually as the wearers muscles get stronger and more familiar with Robohand use.

The Robohand wearer should stretch their wrist bending (flexing) muscles every two hours (holding stretch for 10 seconds) while they are getting used to working the Robohand. Pictures of these stretches are below. This is to keep flexibility in the forearm muscles as they get stronger.*

*If the Robohand wearer's elbow becomes tender, this means the muscle is stronger then flexible and the Robohand should not be used till the pain has stopped. The stretches should be continued till tenderness resolves. Ice pack on elbow is helpful.

(This is the arm being stretched)

